*MANUALE D'USO
USER MANUAL
MODE D'EMPLOI
HANDBUCH
MANUAL DE INSTRUCCIONES*

CE 0123

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IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2002/96/EC:

At the end of its working life, the product must not be disposed of as urban waste. It must be taken a special local authority differentiated waste collection or to a dealer providing this service. Disposing of a household appliance separately avoids possible negative consequences for the environment and health deriving from inappropriate disposal and enables the constituent materials to be recovered to obtain significant saving in energy and resources. As a reminder of the need to dispose of separately the product is marked with a crossed-out wheeled dustbin.

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NEW ASKIR 20 it's a device working with 230V ~ / 50 Hz network electricity, to be used for the nasal aspiration, oral aspiration, tracheal aspiration of the body liquids (mucus or catarrh) in the adult or in the children.
 Easily transportable from one hospital ward to another, particularly suitable for tracheotomized patients, minor surgical applications and post-operative therapy at home.
 Easily portable equipment designed for not continuous use 30 min. ON / 30 min. OFF. Made of highly heat-resistant, electrically insulated plastic material in conformity with the latest European safety standards.
 Supplied with a complete polycarbonate autoclavable jar complete with overflow valve.
 Equipped with aspiration regulator and vacuum indicator located on the front panel.

GENERAL WARNING

READ INSTRUCTION MANUAL CAREFULLY BEFORE USE

ONLY HIGHLY QUALIFIED STAFF USE RESERVED

**THE INSTRUMENT MUST NOT BE DISASSEMBLED
 FOR A TECHNICAL SERVICE ALWAYS CONTACT CA-MI**


IMPORTANT SAFETY RULES

1. On opening the packaging, check the integrity of the appliance, paying particular attention to the presence of damage to the plastic parts, which may make access possible to internal live parts and also to breakage and / or peeling of the power supply cable. **In these cases don't connect the plug to the electric socket. Carry out these controls before each use;**
2. before connecting the appliance always check that the electric data indicated on the data label and the type of plug used, correspond to those of the mains electricity to which it's to be connected;
3. If the plug supplied with the appliance is incompatible with the mains electricity socket, contact qualified staff for replacement of the plug with a suitable type. The use of simple or multiple and / or extension adapters is not generally recommended. Whenever their use is indispensable, use those in compliance with safety regulations, however paying attention not to exceed the maximum power supply limits, which are indicated on the adapters and extensions;
4. Respect the safety regulations indicated for electrical appliances and particularly:
 - Only use original accessories and components;
 - The device can be used only with the bacteriological filter;
 - Never immerge the appliance into water;
 - Position the appliance on flat stable surfaces;
 - Position the device in a way that the air inlets on the back aren't obstructed;
 - Never use the device in environments which have anaesthetic mixtures inflammable with air, oxygen or nitric oxide;
 - Don't touch the device with wet hands and always prevent the appliance coming into contact with liquids;
 - Keep off the reach of children or not capable people without supervision;
 - Don't leave the appliance connected to the power supply socket when not in use;
 - Don't pull the power supply cable to disconnect the plug remove the plug from the mains socket correctly;
 - Preserve and use the medical device in environments protected from atmospheric factors and at a distance from heat sources;
5. For repairs, exclusively contact CA-MI technical service and request the use of original spare parts. Failure to comply with the above can jeopardise the safety of the device;
6. **This medical device must be destined exclusively for the use for which it has been designed and described in this manual.** Any different use must be considered incorrect and therefore dangerous; the manufacturer cannot be considered liable for damage caused by improper, incorrect and / or unreasonable use or if the appliance is used in electrical plants that are not in compliance with the regulations in force;
7. Particular precautions must be made concerning electromagnetic compatibility. The medical device must be installed and used according to information supplied with the accompanying documents;
8. Instrument and accessories discharging must be done following current law regulations in every country of use.

TECHNICAL CHARACTERISTICS

TPOLOGY (MDD 93/42/EEC)	Dispositivo Medico Classe IIa	
MODEL	NEW ASKIR 20	
UNI EN ISO 10079-1	HIGH VACUUM / LOW FLOW	
POWER FEEDING	230V~ / 50Hz	110V~ / 60Hz (no CE 0123)
POWER CONSUMPTION	184 VA	105 VA
FUSE	F 1 x 1.6A 250V	F 1 x 4A 250V
MAXIMUM SUCTION PRESSURE (without jar)	-75kPa (-0.75 Bar) Adjustable from -75kPa (-0.75 bar) to -10kPa (0.10 bar)	
MAXIMUM SUCTION FLOW (without jar)	16 l/min	
WEIGHT	2.5 Kg	
SIZE	350 x 210 x 180mm	
DUTY CYCLE (to 35°C and 110% operating voltage)	30 min. ON / 30 min. OFF	
WORKING CONDITION	Room temperature:	10 ÷ 40°C
	Room humidity percentage:	20 ÷ 85% RH
	Altitude:	0 ÷ 2000m s.l.m.
CONSERVATION CONDITION AND TRASPORT	Room temperature:	-40 ÷ 70°C
	Room humidity percentage:	10 ÷ 95% RH

SYMBOLS

	Class II isolation equipment
CE 0123	CE marking in conformity with EC directive 93/42/EEC
	Warning, consult the instruction manual
	To Preserve in place coolness and dry land
	Conservation temperature: -40 ÷ 70°C
	Type B equipment
	Fuse
~	Alternate Current
Hz	Mains Frequency
I	ON
O	OFF

Guidance and manufacturer's declaration – Electromagnetic Emissions (as request by regulation EN 60601-1-2:2001)		
The surgical aspirator NEW ASKIR 20 is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator NEW ASKIR 20 should assure that it's used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment - guidance
Irradiated / Conducted emissions CISPR11	Group 1	The surgical aspirator NEW ASKIR 20 only used RF energy only for its internal functioning. Therefore its RF emissions are very low and are not cause interference in proximity of any Electronic appliances.
Irradiated / Conducted emissions CISPR11	Class [B]	The surgical aspirator NEW ASKIR 20 can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes.
Harmonic emissions IEC/EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

Guidance and manufacturer's declaration – Electromagnetic Emissions (as request by regulation EN 60601-1-2:2001)		
The surgical aspirator NEW ASKIR 20 is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator NEW ASKIR 20 should assure that it's used in such an environment.		
Immunity Test	Compliance	Electromagnetic environments - guidance
Electrostatic discharge (ESD) IEC/EN 61000-4-2	± 6kV on contact ± 8kV in air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC/EN 61000-4-4	± 2kV power supply	Mains power quality should be that of a typical commercial environment or hospital
Surge IEC/EN 61000-4-5	± 1kV differential mode	Mains power quality should be that of a typical commercial environment or hospital
Loss of voltage, brief voltage interruptions and variations IEC/EN 61000-4-11	5%U _T for 0.5 cycle 40%U _T for 05 cycle 70%U _T for 25 cycle <5%U _T for 5 sec	Mains power quality should be that of a typical commercial environment or hospital If the user of the surgical aspirator NEW ASKIR 20 request that the appliance operates continuously, the use of a continuity unit is recommended.
Magnetic field IEC/EN 61000-4-8	3A/m	The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low.
Conducted Immunity IEC/EN 61000-4-6	3Vrms 150kHz to 80MHz (for appliances that aren't life - supporting)	-
Irradiated Conducted IEC/EN 61000-4-3	3V/m 80MHz to 2.5 GHz (for appliances that aren't life - equipment)	-
Note U _T is the value of the power supply voltage		

ACCESSORIES SUPPLIED

DESCRIPTION
COMPLETE ASPIRATION JAR 1000cc
CONICAL FITTING
TUBES SET 6 mm x 10 mm
ASPIRATION PROBE CH20
ANTIBACTERIAL FILTER

The filter is produced with (PTFE) hydrophobic material which prevents fluids entering the pneumatic circuit. When the filter is wet, it's not possible to use the unit therefore the filter should be changed immediately. In case of possible contamination or discolouration, change the filter immediately. Don't use the suction unit without the protection filter fitted. If the suction unit is used in an emergency or in a patient where the risk of contamination is not known the filter must be changed after each use. Available under request with different versions with complete jar 2000cc. The suction probe must be changed at the end of each application.

CLEANING OF ACCESSORIES

Wearing disposable latex gloves, proceed to clean the accessories as follows:

1. Grip the jar with your hand and turn the cap in a counter-clockwise direction.
2. Remove the red cage covering the float located in the jar cap (conical part with 3 openings)
3. Once the cage is removed, the float with the red seal comes out of the jar cap.
4. **To disinfect the accessories, use either denatured alcohol or hypochlorite-based solutions, easily purchasable from chemist shops.**
5. Once cleaning has been completed, insert the float in the cage, with the seal facing towards the cage opening.
6. Insert the opening of the cage into the seat in the jar cap.
7. To facilitate the cleaning of the jar, fill it 1/3 full of ordinary water. The aspirated liquid, thus diluted, will then be easier to remove.

It is possible to autoclave the accessories cover and jar: insert the parts in autoclave and to effect a cycle of sterilization having care to position turned upside-down the graduated jar (with turned fund toward the tall one). After the sterilization and the cooling to the temperature the components' environment, to verify that this last don't result damaged; disassembled therefore the container for inhaled liquids following the inverse operations to the dismantlement.

The aspiration tubes can be sterilized on autoclave using a sterilization cycle at 120°C.

The conical connector can be sterilized on autoclave using a sterilization cycle at 121°C.

The device is ready for a new employment now.

During the sterilization operations the operator must use disposable gloves in latex, shirts and opportune mask not to enter contact with possible substances contaminating

DO NOT WASH, STERILIZE OR PUT IN AUTOCLAVE THE ANTIBACTERIAL FILTER

MAINTENANCE

The **NEW ASKIR 20** suction equipment does not need maintenance or lubrication.

It is necessary to check functioning and instrument before every use.

Unpack the instrument and **always check** integrity of plastic parts and feeding cable, they might have been damaged during previous use.

Connect cable to electrical network and turn switch on.

Close the aspiration outlet with your finger and with suction regulator in maximum vacuum position check that the vacuum indicators reaches -75 kPa (-0.75 bar) maximum.

Rotate the knob from right to left and check the aspiration regulating control.

The vacuum indicator should go down -25 kPa (-0.25 bar).

Verify that loud noises are not present, these can indicate wrong functioning.

A protection fuse (**F 1x1.6A 250V and F 1 x 4A 250V for voltage 110V/60Hz**) reachable from exterior and it situated in the plug protects the instrument. For use replacing, always check the type and the range indicated.

Fault type	Cause	Solution
1. No aspiration	Jar Cap badly screwed down	Unscrewed the cap, then rescrew it correctly
2. No aspiration	Lid seal not in its seat	Unscrew the cap and insert the seal properly in its seat
3. The float doesn't close	If the cap has been washed, ensure that the float is not partially detached	Insert the float into it's place
4. The float doesn't close	The float it's covered by dirty material	Unscrewed the cap, leave the and put in on autoclave
5. Low suction	Foam inside the jar	Fill the jar to 1/3 full of ordinary water
6. No aspiration due to flow leakage of mucus	Filter blocked	Replace filter
Faults 1 - 2 - 3 - 4 - 5 - 6	None of the remedies has achieved the desired results	Contact the seller or CA-MI After-sales Assistance Service

If the overflow security system it's activated, don't proceede with the liquid aspiration.

If the overflow security system doesn't work there are two cases:

1° case – If the overflow security system doesn't work the aspiration will be stopped by the bacteriological filter who avoid the liquid penetration inside the device.

2° case – If both the security system doesn't work, there is the possibility that liquid comes inside the device, in this case return the device to CA-MI technical service.

BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING, PLEASE CONTACT CA-MI TECHNICAL SERVICE.

CA-MI DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE CHECKING, APPEARS TO BE TAMPERED.

INSTRUCTIONS

- Connect the short silicon tube **1**, with antibacterial filter **8**, to the suction connector **2**. The other tube, with one end connected to the filter must be connected with the other end of jar's lid connector **3** where has been fixed the red float (security float). When the 90% of the volume of the jar is reached there is the activation of the security float (the float close the aspiration connector on the jar) to avoid liquid penetration inside the device. The device must be used on a plan of horizontal operation.
- Connect the long silicon tube **4** to the other jar's lid connector **5**.
- Connect the other end of the long silicon tube **4** to the probe plastic connector **6** then connect the suction probe to it.
- Unscrew the jar's lid and fill the jar 1/3 full of ordinary water (this for an easy cleaning operations and an rapid reaching of the functionality vacuum) then rescrew the lid on the jar correctly.
- Connect the power cord to the device then connect the plug to the electrical mains supply.
- Push switch **7** on position **I** to start suction
- Once finished push switch on **O** position and unplug.

RULES FOR RETURNING AND REPAIRING

COMPLYING WITH THE NEW EUROPEAN RULES, CA-MI INDICATES THE IMPORTANT POINTS TO PROTECT INSTRUMENT AND OPERATORS HYGIENE. THESE RULES MUST BE RESPECTED IN ORDER TO GUARANTEE HYGIENE AND SAFETY TO ALL THE PEOPLE OPERATING WITH THE INSTRUMENT TO OBTAIN QUALITY AND WELL BEING.

Every returned instrument will be hygienically checked before repairing. If CA-MI finds instrument not suitable for repairing due to clear signs of internal or external contamination, the same will be returned to customer with specification of NOT REPAIRED INSTRUMENT, accompanied by an explanation letter.

CA-MI will decide if contamination is due to bad functioning or misuse. If contamination is due to bad functioning, CA-MI will substitute the instrument, only if a SALE RECEIPT and STAMPED GUARANTEE accompany the same.

CA-MI is not responsible for contaminated accessories, they will be substitute at customer's expenses.

For this reason it is **COMPULSORY** to carefully disinfect the external part of the instrument and accessories with a cloth soaked in methylated spirits or hypochlorite-based solutions. Put the instrument and accessories in a bag with indication of disinfecting. We also request to specify the kind of fault, in order to speed up repairing procedures.

To this end, please read the instructions carefully in order to avoid damaging the equipment through improper use.

Always specify the fault encountered so that CA-MI can establish whether it falls into the category of the faults covered by the guarantee.

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