


Operation manual Service manual

VETLux 150

150 Watt Halogen light source with connector
and air pump for flexible endoscopy



Eickemeyer - Medizintechnik für Tierärzte

Eltastr. 8 - 78532 Tuttlingen, GERMANY
Tel. +49 (0) 7461 / 96580-0 Fax: +49 (0) 7461 / 96580-90
Internet: www.eickemeyer.de

This manual contains information protected by copyright laws. All rights reserved. Manual may not be copied or duplicated, in its entirety or in part. We reserve the right to make technical changes on device prior notice.

Copyright Eickemeyer
10.2005

SERVICE ADDRESS:

Eickemeyer Medizintechnik für Tierärzte
Eltastr. 8
D-78532 Tuttlingen
GERMANY

Telephone: +49 (0) 7461 / 96 580-0
Fax: +49 (0) 7461 / 96 580-90



Table of contents

GENERAL ADVISES	4
DATA OF THE EQUIPMENT	4
RIGHTS	4
WARRANTY	4
SERVICE, REPAIRS AND MODIFICATIONS	5
LIABILITY	5
DISPOSAL.....	5
SAFETY REFERENCE / PLACE THE EQUIPMENT	6
NORMAL USE	6
USER QUALIFICATION	6
UNPACKING / ITEMS INCLUDED	6
SAFETY PRECAUTIONS OF INSTALLATION.....	6
STORAGE AND OPERATION CONDITIONS.....	6
SIGNS & SYMBOLS	7
SYMBOLS	7
SIGNS.....	7
DESCRIPTION OF THE EQUIPMENT	8
OPERATING ELEMENTS	9
FRONTPANEL	9
REAR PANEL	9
CONNECTING AND OPERATING	11
CONNECTING THE EQUIPMENT	11
OPERATING THE EQUIPMENT.....	12
SERVICE MANUAL	13
GENERAL MAINTENANCE AND REPAIR ADVICE	13
SERVICE INTERVALS	13
EXCHANGE OF THE MAINS FUSES	13
EXCHANGE OF LAMP	14
CLEANING / DISINFECTING.....	14
WIRING DIAGRAM	15
TECHNICAL DATA	16
SPARE PARTS	16
TABLE ‘TECHNICAL SERVICE-INFORMATION’	17
APPENDIX	18
ELECTROMAGNETIC COMPATIBILITY (EMC).....	18

General advises

Like all of our other products, this product is the result of years of experience and great care in engineering and manufacture.

This manual is destined to learn you understanding the function and the operation of your equipment.

Before you switch on the equipment for the first time, please thoroughly read this manual and pay special attention to all safety instructions, so that endangering for the user and the patient is precluded.

Please always store this manual with the equipment.

Before use, read this manual thoroughly.

An insufficient understanding of the dangers, warnings, cautions, and informations in this manual can result in death, serious injury, or equipment damage.

This product complies with the requirements of Directive 93/42/EEC concerning medical devices.



Data of the equipment

The type label (rear of unit) contains technical data, type and serial number of your unit. Please always indicate these data when ordering spare parts or in case of any question.

Please enter here the technical data of your device! →

- Serial No.: _____
- Type: _____
- Date: _____
- Class: _____
- Hz: _____
- Amp.: _____
- Volt: _____

Rights

All rights on this user manual, especially the rights of duplication and publication and the rights on translation are reserved. No part of this user manual may be reproduced by any means (by photostatic copy, microfilm or other methods) without preceding written consent of the manufacturer or be reproduced, multiplied or published by means of electronic data processing. The information given in this user manual can be changed or extended without notice and do not represent any liability. Errors and technical changes excepted.



Warranty

1 year according to our warranty conditions.

Opening the equipment or performance of any repairs or modifications of the equipment by unauthorised persons shall relieve the manufacturer of any liability for its performance. Any such opening, repair or modification performed during the warranty period shall void all warranty.

Wear parts are not included in the warranty.

The firm of the manufacturer shall be liable for failure or deterioration in the safe operation, operational reliability and performance of this equipment only subject to the conditions, that all assembly operations, system expansions, readjustments or repairs to same have been performed by a person or persons duly authorised by the manufacturer, that all electrical installations at the location of us meet applicable national and local electrical codes and that the instrument has been used in accordance with its operating instructions at all times.

Service, Repairs and Modifications

In conformity with the international safety regulations valid for medical devices, all activities such as check ups, repairs, modifications, calibrations etc. may only be carried out by the manufacturer or by explicitly authorized personnel.

All services carried out must be entered in the „Technical Service Notes“ at the back of the user manual.

Liability

As manufacturer of the device, we only consider ourselves liable for safety, reliability and performance of the unit, if

- assembly, re-adjustment, modifications or repairs are performed by persons authorised by us
 - the electric installation of the respective room corresponds to the regulations of VDE 0107
- the instructions found in the user manual are strictly observed when operating the unit.

Disposal

In accordance with European Directive 2002/96/EC on Waste Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately.

Refer to your distributor resp. the manufacturer for return and/or collection system available in your country.



Safety reference / Place the equipment



Normal use

The equipment may only be used with accessories, wearing parts and disposable items, which have been designated by the manufacturer suitable for the instrument or the safety use of which is proven.



User qualification

The equipment may only be used by persons, who have a corresponding specialised qualification and who have been instructed in use of the equipment.

It is the user's responsibility to make sure, the equipment is safe and operates properly before using the equipment.

Unpacking / items included

Carefully unpack the equipment and accessories and remove it from their packing.

Check for missing items and evidence of shipping damage.

File any complaints with the manufacture or supplier immediately.

Retain the original packing materials for later use. These can come in handy, when the equipment must be transported.

Please verify immediately after having unpacked the equipment, whether the delivery is complete. The standard extent of delivery includes of the following:

- Control unit
- Power supply cord
- Manual



Safety precautions of installation

Always place the equipment on a solid base.

The equipment may be used only in rooms having electrical installations conforming to applicable national, state and local electrical codes.

The unit must be joined to the central potential equalisation of the operating theatre or of the equipment trolley by means of a grounding cable.

The device must be connect to line voltage using the delivered protectively earthed power supply cord.

Storage and operation conditions

Storage temperature: -20°C to +60°C

Operation temperature: +10°C to +40°C

rel. humidity: storage: 10% to 90%
 operation: 30% to 75%

Air pressure: storage: 600 mbar to 1300 mbar
 operation: 700 mbar to 1060 mbar

Signs & symbols

Symbols

Attention, important note!



Safety note



Service



Signs

Please read the enclosed instructions !



Unit model BF



Beware of dangerous electrical voltage !



Connection for ground potential



Alternating voltage



Description of the equipment

The HL 150 is a powerful 150W-Halogen lightsource for all applications with fibre-optic illumination. The lightintensity makes the lightsource usable for almost all endoscopic disciplin. The output intensity can be easily controlled via control knob without influence on the quality of light.

The built in soft-start-electronic protects the lamp against high currents in the switch on/off moments and extends the lamp-lifetime.


The lightsource comes equipped with a built-in insuflation pump for use with gastroscopes and other flexible endoscopes. The pump can be switched on/off on the frontpanel.

There are adapters available for gastroscopes made by either Olympus, Fujinon or Pentax. The pump can be separately switched on and off.

The equipment is desinged by the latest standards of security for medical equipment and complies with the demands of CE.


Operating elements

Frontpanel



- (1) Mains switch
- (2) Switch Air pump
- (3) Intensity adjustment
- (4) Light guide connector

Rear panel



- (5) Connector potential ground
- (6) Line connector
- (7) Fuse holder

Mains switch

The control unit is turned on by switching the mains switch.

The mains switch has two different switching positions:

- I switched on
- O switched off

When the control unit is switched on, this is indicated by the up-light green lamp inside the switch.

Intensity control knob

Use this knob to control the output intensity. Turning the knob to the left decreases intensity, turning to the right increases intensity.

Pump switch

The integrated pump is switched on and off by using the corresponding switch.

Switch on the pump only right before the use of it and switch it off right after.

The pump is maintenance-free. Make sure, no dirt and liquids can get into the adaptor.

Connector for fibrescope

For lightguide connection, there are interchangeable adapters available for all brands of lightguides. For gastroscopic applications, the lightsource comes equipped with a special adapter with integrated air supply for your gastroscop brand.

Using the wrong lightguide or adaptor can damage the lightsource and the accessories and threatens the patient.

Potential ground connector

Basically, the control unit is protectively earthed by the 3-pin power supply cord when it is connect to a protectively earthed wall socket, as prescribed.

When running the equipment in rooms which comply to class 1 or 2E according to MedGV, the control unit must be joined to the central potential equalisation of the operating theatre or of the equipment trolley by means of a grounding cable.

Fuse holder

It contains the mains fuses. You have to control whether your mains voltage corresponds with the selection shown in the window.

Line connector

The plug of the power supply cord is connect to the mains terminal device. Use only the delivered supply cord.

Connecting and operating

Connecting the equipment

Before you connect the mains plug, check on the back of the equipment that the voltage indicated (in the square panel above the mains socket) is the correct one.

voltage = 230VAC → indicator '230'

voltage = 115VAC → indicator '115'

If the incorrect voltage is indicated then the equipment may under no circumstances be connect. Before connecting a light cable into a light source adapter of the equipment, please ensure that you have the correct plug and adapter. A plug that is too long or too thin can for example be pushed too far into the equipment, and can damage the sensitive diaphragm or lens inside. This will result moreover in considerable loss of light, in the same way as a plug that is too short, or a plug of the wrong diameter, because of the wrong position of the contacts at the light entry.



All connection have to be done before you switch on the unit !

Connect the power supply cord !

Use the delivered protectively earthed power supply cord to connect the control unit to the mains. Make sure that on comes indicated mains voltage is correct. Attach the equipment only to agrounded protective contacz socket.

Connect a fibrescope !

Verify, whether the installed adapter matched the fibrescope you want to connect. Plug in the light plug of the fibrescope. Thus the air connection is automatically fixed. Connect the water bottle to the fibrescope (if necessary). There is a holder for the water bottle installed at the left side of the housing.

Connect the potential equalisation conductor !

Join the terminal device for potential equalisation on the rearpanel with the central potential equalisation of the operating theatre or of the equipment trolley.

Operating the equipment

Before connecting a light cable into a light source adapter of the equipment, please ensure that you have the correct plug and adapter. A plug that is too long or too thin can for example be pushed too far into the equipment, and can damage the sensitive diaphragm or lens inside. This will result moreover in considerable loss of light, in the same way as a plug that is too short, or a plug of the wrong diameter, because of the wrong position of the contacts at the light entry.

After switching on a light source and the ensuing ignition of the lamp, the equipment should remain switched on for at least $\frac{1}{4}$ hour. A shorter shining period will considerably shorten the life expectancy of the lamp! After switching off a light source however, it may be immediately switched on again. A waiting period or cooling off period is not necessary.

After connecting the mains, the grounding conductor, the light cable and the endoscope, the equipment is ready for use, and can be switched on at the mains switch. The green lamp inside the switch lights up. Now you can begin working with the lightsource. Set the desired intensity by turning the corresponding knob. There are no other control elements available to use during operation.

Set the pump to on or off by switching the corresponding switch. Please switch on the pump only when necessary

Service manual

General maintenance and repair advice



The instructions and information given in this chapter are only for instructed personnel, who are aware of the safety precautions necessary for repair and maintenance of medical electronic devices.

The manufacturer refuse any liability for unauthorised repair and modification.

The manufacturer will provide those circuit diagrams, itemised parts listings, descriptions, sets of adjustment instructions and other items of available documentation to suitably qualified user personnel duly authorised by the manufacturer for their use in repairing those components of the equipment that have been designated by their respective manufactures as repairable.

Only the supply of such technical documentation relating to the equipment shall not be construed as constituting manufacturer's authorisation of user's personnel, regardless of their levels of technical training, to open or repair the equipment.

Explicitly exempted here from are those maintenance and repair operations described in this manual.

Service Intervals

The equipment should be controlled on save function by authorized service personal at least one time the year. These controls should be done every time too, when a lamp is exchanged.

Exchange of the mains fuses



The mains fuses are located on the rearpanel of the control unit, right above the mains terminal device in a small drawer. If you need to exchange the mains fuses, proceed as follows:

- **PULL OUT THE MAINS PLUG!**
- Loosen the drawer by unfastening the two clamps located to the left and to the right of the drawer with a peaked tool and pull out the drawer.
- Take out the fuses.
- Check the fuses. A blown fuse is indicated by the blackened glass cylinder or the visibly melted fuse conductor. If necessary, check the fuse with an ohmmeter.
- Install the corresponding fuses.
- Re-install the fuse-drawer.

Switch on the equipment again. If you have exchanged a defective fuse against a new one and the fuse blows again, the unit has an error. In this case, you must return the device to your dealer for testing and repair.



Exchange of lamp


The only part subject to wastage in the equipment are the lamp. The manufacturer of the lamp estimates the minimum life expectancy as 50 hours. After that, the quoted light wattage reduces drastically and there may be difficulties in igniting. In order to secure safe operation of the equipment, the lamp should be changed after expire of this operation time.

The lamp must be changed for a new one as following:

- ***PULL OUT MAINS PLUG!***
- Turn the light source button up!
Loosen the screw on the bottom and remove the lamp cover.
- Attention: If the machine has been operating shortly beforehand the lamp could be very hot.
- All remaining current in the electronics will be, in a few seconds after switching off the equipment, completely safely discharged, so there is no danger in handling the equipment here.
- Pull out the lamp upwards out of its holder.
- Hold fast the reflector of the lamp and remove the cable from the lamp.
- Install the cable now to the new lamp.
Insert a new lamp into the holder and push it downwards.
- CHECK, the little nose of the lamp is in the right pocket.
- Re-install the cover and tighten it with the screw.
- Observe the grounding cable!

Further servicing, which should be done in conjunction with lamp replacement

Dust should be removed from the equipment with a vacuum cleaner using a small nozzle attachment. Dust can sometimes cling quite firmly to the ventilator blades. This should then be cleaned with a cloth and a little alcohol/spirits. Likewise dust settles on the heat protective filter (with lens, optional), this should also be cleaned with a soft cloth or blotting paper and pure alcohol or spirits.




Cleaning / Disinfecting

NOTE: PULL OUT THE MAINS PLUG!

All parts of the outer surfaces of the equipment are totally insensitive to all the usual cleaning and disinfecting materials, so that you can use any of these without limitation. Apply liquids using a soft cloth or soft blotting paper, in order to avoid scratches on the surfaces and in order to be able to control the amount of liquid.

With flammable liquids like alcohol especially, you should apply with a cloth. Do not let any liquid get into the equipment. After cleaning with flammable liquids, leave the equipment to dry for one hour, before it is switched on again. There is danger for example that an alcohol-air explosive mixture could form after cleaning.

Wiring diagram



TITLE: HL 150 mini Wiring Diagram	
Document Number:	REV:
Date: 29.10.2003 08:15:12	Sheet: 1/1

Technical Data

Power supply	230 VAC \pm 10% or 115 VAC \pm 10%
Lamp	Halogen reflector lamp, 15V, 150W
Power consumption	160 W
Fuses	2 fuses 5x20 mm 230 VAC: 1,2,5 A 115VAC: 2,5A
Dimensions	245 x 120 x 200 mm (W x H x D)
Weight	3,5 kg
Class / Type	1 / BF
Manufactured and tested acc.	IEC 601-1 / CE; EN 60601-1; 93/42/EEC

Spare parts

Lamp	Halogen reflector lamp 15 V, 150 W OSRAM HLX 64634, EFR
Mains fuses:	fine fuses, 5x20mm, delay acting for 230VAC: 1,25A for 115 VAC: 2,5mA

Appendix

Electromagnetic Compatibility (EMC)

Precautionary measures

Electromedical devices are subject to special precautionary measures concerning electromagnetic compatibility (EMC). This device is to be used for the purposes described in the operation manual and has to be installed, set up, and operated in compliance with the EMC guidelines.

Impact of mobile and portable RF communication devices

The emission of high frequency by mobile communication devices may impact the function of the electromedical device. Operating such mobile communication devices (e.g. cell phones, GSM phones) in the proximity of the electromedical device is prohibited.

Electrical connections

Connections between such plugs and sockets may not be established without first implementing ESD precautionary measures.

ESD precautionary measures

The following are ESD precautionary measures:

- Connect all electrical equipment to be connected to the device to a potential equalisation system (via PE).
- Use only the equipment and accessories mentioned in the operation manual.

The staff have to be informed about and trained in ESD precautionary measures.

The VETLUX 150 should be operate in an enviroment as stadet below. The user have to make sure that the VETLUX 150 will be operate in such enviroment.

Manufacturer declaration – Electromagnetic application

Noise emmission measurement	Correspondence	Electromagnetic enviroment - Guide
HF- Emission acc. to CISPR 11	Group 1	The VETLUX 150 used HF-energy solely for internal functions. Therefore, the HF-emission is very low and interferences with adjacent electronical equipment are improbable.
HF- Emission acc. to CISPR	Class B	The VETLUX 150 is adapted for use in all facilities including living areas and such facilities with public utility provider.
Harmonics acc. to IEC 61000-3-2	Class A	
Voltage fluctuations / Flicker acc. to IEC 61000-3-2	unapplicable	


Manufacturer declaration – Electromagnetic noise immunity

Noise immunity tests	IEC 60601 – test level	Correspondence level	Electromagnetic enviroment - Guide
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be consist of concrete or wood or furnished with ceramic tiles. If the floor furnished with synthetic material, the relative air humidity should be 30 % average.
Quick transient electrical disturbances / Bursts acc. to 61000-4-4	± 2 kV for mains lines ± 1 kV for Input and output lines	± 2 kV for mains lines ± 1 kV for Input and output lines	The quality of mains voltage should be comply typical mains voltage of business or hospital enviroment.

Surge voltage acc. to IEC 61000-4-5	± 1 kV push-pull voltage ± 2 kV push-push voltage	± 1 kV push-pull voltage ± 2 kV push-push voltage	The quality of mains voltage should be comply typical mains voltage of business or hospital enviroment.
Acc. to IEC 61000-4-11	< 5% U _T (> 95% U _T collaps)for ½ period < 40% U _T (> 60% U _T collaps)for 5 periods < 70% U _T (> 30% U _T collaps)for 25 periods < 5% U _T (> 95% U _T collaps)for 5 sec.	< 5% U _T (> 95% U _T collaps)for ½ period < 40% U _T (> 60% U _T collaps)for 5 periods < 70% U _T (> 30% U _T collaps)for 25 periods < 5% U _T (> 95% U _T collaps)for 5 sec.	The quality of mains voltage should be comply typical mains voltage of business or hospital enviroment. If the user of the VETLUX 150 demand continued functions also by appearance of interrupts, it will be recommended to operate the VETLUX 150 by non-interruptable power supply.
Magnetic field by supply frequency (50 / 60 Hz) acc. to IEC 61000-4-8	3 A / m	3 A / m	Magnetic fields by line frequency should be correspond to typical values you will find in business and hospital enviroments.
Note: U _T is the line AC voltage before using the test level.			

Manufacturer declaration – Electromagnetic noise immunity

Noise immunity test	IEC 60601 – test level	Correspondence level	Electromagnetic enviroment - Guide
			Portable and mobile radio transmitter should not be use within the recommended safety distance to the VETLUX 150 (including the lines)! The recommended safety distance will be calculate by the appropriate transmitting frequency formula. Recommended safety distance:
Guided HF-disturbances acc. to IEC 61000-4-3	3 V _{eff} 150 KHz to 80 KHz	3 V _{eff}	d=1,2 √P
Radiated HF-disturbances acc. to IEC 61000-4-3	3 V / m 80 MHz to 2,5 GHz	3 V / m	d=1,2 √P 80 MHz to 2,5 GHz
			d=2,3 √P 80 MHz to 2,5 GHz
			P – transmitter wattage rating (W) acc. to manufacturer (transmitter)data d – recommended safety disdnce in meters (m) The electrical field strength of fixed radio transmitter is lower than the correspondence level ^b according to a

			fieldwork ^a . In the environment of equipment which bore the following pictograph, it is possible to appear interferences. 
Note 1: At 80 MHz and at 800 MHz the higher result is valid.			
Note 2: This guideline may not be the case for all situations. The spread of electromagnetic waves will be affected by absorption and reflexion of buildings, objects and humans.			
a	The electric field strength of fixed transmitter, like wireless telephones, mobile radio services, AM- and FM broadcast transmitter and TV transmitter, couldn't be predetermined correctly. In order to ascertain the electromagnetic environment at fixed HF-transmitter, it is commendable to inspect the location. If the ascertained field strength at the VETLUX 150 location exceeds the specified correspondence level, it's may be necessary to adopt other means, e.g. reorientation or removal the VETLUX 150.		
b	By means of the frequency range from 150 kHz to 80 MHz the field strength is less than 3 V/m		

Commendable safety distances between portable and mobile communication appliances and the VETLUX 150

The VETLUX 150 is intended for operation in electromagnetic environments, where radiated HF-disturbances will be controlled. The customer resp. the user of the VETLUX 150 can help to prevent electromagnetic disturbances by compliance minimum distances between portable and mobile HF communication appliances (transmitter) and the VETLUX 150.

Commendable safety distances acc. to the maximum power output of the communication appliance:

Wattage rating of the transmitter (W)	Safety distance acc. to transmitting frequency (m)		
	150 kHz to 80 MHz $d=1,2 \sqrt{P}$	80 MHz to 800 MHz $d=1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d=2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitter, whose wattage rating isn't stated above, the safety distance can be calculated by using the equations in the respective column.

P is the wattage rating of the transmitter acc. to the specifications of the transmitter manufacturer.

Note 1: To calculate the commended safety distance of transmitter, with a frequency range from 80 MHz to 2,5 GHz, was an additional coefficient of 10/3 used. The additional coefficient should reduce the probability of disturbances by portable / mobile communication appliance, which are unintentional brought in the patient range.

Note 2: This guideline may not be the case for all situations. The spread of electromagnetic waves will be affect by absorption and reflexion of buildings, objects and humans.