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SONOREX DIGITEC Serie F

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Operating Instructions

SONOREX DIGITEC Serie F

ultrasonic units for sample preparation

Units, accessories and agents have to be strictly used according to the instruction manual or the corresponding product information.

Keep units and agents away from children and from persons not instructed for correct use by means of these instructions.

In case of improper use, the manufacturer does not take any responsibility for safety and proper functioning.

Unauthorized changes/modifications will result in expiration of the CE-conformity.

If service is required, please contact your dealer or the address mentioned below.

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BANDELIN electronic GmbH & Co. KG

Copyright & limit of liability

This operating instruction may not be reproduced, either in full or in extracts, without prior approval of BANDELIN electronic GmbH & Co. KG, hereinafter referred to as BANDELIN.

BANDELIN accepts no responsibility and liability for damage caused by improper handling or usage contrary to the intended purpose.

The German original is the binding version of this document. Any difference in the translation is not binding and has no legal effect. In case of any discrepancies between the translation and the original, the original version has priority.

The documentation was created with great care. Liability for indirect and direct damage arising because of incomplete or erroneous information in this documentation as well as in delivery and usage is excluded.

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Contents

1	Strictly observe the following general instructions for safety and correct operation!	14
2	Zweckbestimmung	15
3	Supply	15
3.1	Available Accessories	15
4	Start	16
4.1	Degassing	16
4.2	Probe preparation	16
4.3	Explanation of operating elements and start-up	17
4.4	Ultrasound operation	18
4.5	Leads for operation and special funktions	19
4.6	Test routine for operating functions	19
4.7	Emptying the ultrasonic tank	20
5	Cleaning and care	20
6	Technical data	21
7	Disfunctions	22
	Notes	23

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Strictly observe the following general instructions for safety and correct operation!

- Keep units and agents away from children and also away from persons not instructed for correct use by means of these instructions.
- Do not immerse parts of the body (hand, feet) nor creatures or plants into the ultrasonic tank. During ultrasonic operation do not dip your hands into the liquid.
- Avoid any contact of your eyes and/or your skin with the agents.
- Do not ingest nor breathe in the agents.
- Never run the unit without liquid!
- Do not fill hot water into the ultrasonic tank. The filling temperature must not exceed 50 °C.
- The contact liquids in the ultrasonic oscillating tank may not exceed a maximum operating temperature of 100 °C.
- Do not put parts into the tank with aggressive soiling such as acids, chloride ions or similar substances!
- Do not fill the stainless steel oscillating tanks with aggressive agents such as acids or chemicals which contain or release chloride ions (some disinfectants, dish-washing detergents, household cleaners, saline solutions).
- Distilled or deionized water without additives is unsuitable for ultrasonic operation directly in the oscillating tank. They lack the necessary cavitation properties to form a uniform ultrasonic field. Irregular cavitation can damage the oscillating tank or the objects being cleaned. (You may use distilled or deionized water in the inset beakers and non-perforated insert tubs.)
- Do not fill the stainless steel oscillating tanks with combustible liquids.
- Use aggressive cleaning agents only in inset beakers or non-perforated plastic insert tanks.
- After removing the packing, check the unit's casing and control elements for any signs of damage. In such case, do not connect the unit to the electrical power supply, but report the damage immediately and in writing to the carrier and your supplier. Save the original packing.
- Place the unit in a horizontal position on a dry, solid base.
- Plug the unit only to a socket with grounded safety contact.
- Keep the unit's surface and control elements clean and dry!
- Replace defective parts only with original **SONOREX** parts!
- Ultrasonic operation produces cavitation noise. If working within distances of 2 meters from the unit for longer periods of time, be sure to wear ear plugs to avoid damage to your hearing!
- Do not place anything on the bottom of the oscillating tank. Use appropriate **SONOREX** accessories.
- The ultrasonic energy in the oscillating tank heats up the liquid! Check the temperature regularly when sonicate temperature-sensitive probes.
- Do not move or carry units filled with liquid.

The manufacturer does not take any responsibility for damage of equipment caused by inappropriate chemicals.

Principle of ultrasonic production

PZT oscillators underneath the tank convert electrical energy into mechanical vibrations.

SONOREX units produce oscillations of 35 kHz in the liquid. Tiniest vacuum bubbles are formed which then implode (cavitation).

2 Purpose

In-vitro diagnostic lab process (sample preparation)

The ultrasonic unit can be used for in-vitro diagnostic sample preparation or laboratory processes, where important medical information can be received by sonication of organic as well as inorganic material. In this case these ultrasonic baths must be treated as medical devices complying with the directive for in-vitro diagnostics 98/79/EG.

In case of improper and irregular disinfection and cleaning a micro-biological contamination caused by settlement of microorganism can occur, especially at the rim of the tank or in the drain section. Cross-contamination might be a result. Therefore, the surfaces of the tank and the housing must regularly be disinfected and cleaned.

Disinfection and cleaning must be carried out regularly by the user in compliance to the hygienic plan with surface disinfectant agent certified by VAH or another efficiently classified agent.

3 Supply

Ultrasonic unit DT 100 F,	Flask holder GL100 F,	250 ml TICKOPUR R 33
Ultrasonic unit DT 510 F,	Flask holder GL 510 F,	250 ml TICKOPUR R 33
Ultrasonic unit DT 1028 F,	2 Flask holder GL 510 F,	250 ml TICKOPUR R 33

3.1 Available accessories

Accessories	Code No.	DT 100 F	DT 510 F	DT 1028 F
EK 10	051	●	●	●
EK 25	053	●	●	●
EK 50	055	●	●	●
EK 100	057	●	●	●
GL 100 F	3261	●	-	-
GL 510 F	3262	-	●	● (2 x)
D 3 P	3214	●	-	-
D 10 P	3246	-	●	-
D 28 P	3258	-	-	●
TICKOPUR R 33	830 / 1 Litre 883 / 2 Litres 831 / 5 Litres	●	●	●

Clamp for laboratory EK 25 to fix laboratory flasks to the flask holder up to 100 ml

Flask holder GL 100 F stainless steel, to carry the clamps EK with flasks

Lid D 3 P plastic, to cover the oscillating tank (without flask holder)

4 Start

- Make sure that the ball valve of the unit is closed!
- Fill the oscillating tank up to the filling mark.
- Connect the unit to a grounded socket.

4.1 Degassing

Degassing the liquid enhances the transmission of ultrasound - see also lead for DEGAS function under chapter 4.5.

Contact liquids contain dissolved gases (for example, oxygen). Sonicate cleaning and contact liquid, which have just been poured into the oscillating tank or have been left in it for a longer period of time, for approx. 5 to 15 minutes before starting the cleaning process. During the degassing process the cavitation noise changes. As the loud degassing noise decreases towards the end of the process, the unit appears to run more quietly. Lower noise levels indicate the end of the degassing process and an improved transmission of ultrasound. They do not mean a reduction of ultrasonic energy.

4.2 Probe preparation

INDIRECTLY in laboratory flasks

⇒ **SONOREX** Accessories examples


- Fill the oscillating tank with a tenside-containing contact liquid (for example TICKOPUR R 33 diluted with water) to achieve a good transmission of ultrasound.

Dosage of the contact liquid:

	Filling capacity	TICKOPUR R 33 1%		
DT 100 F	1,0 Litre	990 ml	+	10 ml
DT 510 F	2,5 Litre	2,48 l	+	25 ml
DT 1028 F	5,8 Litre	5,74 l	+	58 ml
		Water	+	agent

When using very calciferous tap water for preparation of the contact liquid, turbidity or precipitation are possible. To prevent it, deionized water must be used.


- Laboratory flasks can be fixed with spring clamps EK to the flask holder GL. First screw the spring clamps tightly on the flask holder.
- Laboratory flasks must be immersed in the contact liquid to a depth of at least 2 cm. See to the tank filling to 2/3.
- Several laboratory flasks with different samples can be sonicated at the same time.

Note:

A direct cleaning in the flat devices (DT ... F) is **not permissible**.

4.3 Explanation of operating elements and start-up

symbol	operation
	Switch the unit on or off (standby). Time scale shows last value. Actual temperature is indicated through blinking LED.
	Time setting through keypress. Permanent operation= ∞.
	Start and Stopp of ultrasound through keypress.

4.4 Ultrasound operation

- Switch on the unit (press key On/Off) – yellow LED of last selected time on the time scale lights up.

Start without modification of time:

- By pressing the key Start/Stop, ultrasound is being switched on.
This is indicated on the time scale by means of a running light of the balance time.
- By pressing the key Time, ultrasound output can be interrupted (pause).
- By pressing the key Time again, ultrasound output will be continued (balance time).

Start with modification of time:

- Select time - by pressing the key Time, the duration will be gradually increased.
Numeral (1 ... 30) = lighting, yellow LED
Continuous operation (∞) = lighting, green LED
- By pressing the key Start/Stop, ultrasound is being switched on.
This is indicated on the time scale by means of a running light of the balance time.
- By pressing the key Time, ultrasound output can be interrupted (pause).
- By pressing the key Time again, ultrasound output will be continued (balance time).
- After run off of selected time, the unit stops and LED of previously selected time lights up.
- Switch off the unit by pressing the key On/Off.
The time last selected will be stored.

Attention! Observe leads for operation and special functions – chapter 4.5

4.5 Leads for operation and special functions

Time scale / duration:

- LED of selected time is lighting constantly, in case of continuous operation (∞) the top (green) LED lights up.
- Ultrasound activity is indicated through LED blinking in second cycle (running light downwards). Running light always starts at the LED for remaining time.

Time scale / DEGAS:

- For degassing (DEGAS) before operation, switch on the DEGAS function. If necessary, select required duration with the key Time. Then press the key Start/Stop for at least 2 sec. Early switch-off through re-pressing the Start/Stop-key.
During degassing, the top green LED is blinking in addition to the time lapse (└┐└┐)
- Switch between ultrasound and degas: if the Start/Stop-key is pressed a long time with ultrasound on, the ultrasound is first turned off and then reactivated with the degas function after about 2 sec.

Time scale / continuous operation:

- In order to avoid inadvertent switching on of continuous operation, continuous operation can be deactivated:
 - Pull the mains plug.
 - Press the Time key and simultaneously pocket the mains plug.
The yellow LED 1 min lights up for confirmation.The function is again activated in the same way. The green LED continuous operation (∞) lights up for confirmation.

Attention:

For safety reasons the unit will be automatically switched off, if no key has been pressed longer than 12 hours.

Leads regarding temperature

- Warmed-up liquids improve the transmission of ultrasound.
- Ultrasonic energy heats up the liquid. Continuous irradiation and/or the covering of the oscillating tank may cause the temperature of the liquid to rise more quickly.
- When reaching a critical temperature, the ultrasonic power is reduced to protect the electronic components inside the ultrasonic device. This way, the internal temperature is prevented from further rising.

4.6 Test routine for operating functions

A test routine can be started in order to carry out an internal function test:

For this, the unit must be switched off. With the key Start/Stop pressed switch on the unit with the key On/Off.

One after the other, all LEDs will light up for one third second. Afterwards, the last selected values will reappear.

With this, the test has been completed successfully.

In case of deviations, the unit has to be sent in for inspection/repair.

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
Do not place the unit into a sink.


Units with drain

- Ball valve
Turn the lever to the outflow direction.

Units without drain

Pour out the liquid over the rear left tank corner. Avoid splashing from below. Pay attention to the cable connection.

5 Cleaning and care

- Be sure to pull the plug before cleaning the unit.
- Do not shower the units or immerse them in water; instead, wipe them with a damp cloth, in- and outside.
- Rinse out the stainless steel tank frequently and wipe it dry.
- Do not use steel wool, scrapers or graters.
- If rings should appear in the stainless steel tank after a longer period of operation, remove them with a customary stainless steel cleaner without scouring additives.
- Metal parts remaining on the stainless steel surface, as well as particles of rust, can penetrate the passive protective layer of the stainless steel, thereby “activating” it and causing it to rust. This rust provokes pittings on the stainless steel surface. Therefore, remove all remaining metal parts such as screws, metal chips and the like; remove minor rust spots immediately with a soft cloth and customary stainless steel cleaners without scouring additives.

SONOREX units do not require maintenance.

6 Technical data

SONOREX compact units are RFI-proofed and C € marked.

Compliance with the limits according to EN 55011, EN 61000-6-1/08.2002 and EN 60601-1-2.

Nominal voltage: 230V~ 50/60 Hz, (115 V on request)
cable length 2 m

High frequency: 35 kHz

Degree of protection: IP 33

Oscillating tank: stainless steel

Type	Code No.	Oscillating tank interior (l x w x d)	Filling capacity	Drain ball valve	ultrasonic peak output*	HF-output	Current consumption	Generator fuses
		mm	l		W	W _{eff}	A	
DT 100 F	3241	240 x 140 x 65	1,0	-	240	60	0,4	T1A
DT 510 F	3242	300 x 240 x 65	2,5	G ½	560	140	0,7	F3,15A
DT 1028 F	3243	500 x 300 x 65	5,8	G ½	1280	320	1,4	F4A

* To achieve an improved efficiency the ultrasound is modulated whereby in combination with SweepTec a four times times value of the HF-output are received as ultrasonic peak output.

Environmental conditions in accordance with EN 61 010-1 (IEC 1010-1)

Degree of soiling: 2 according to IEC 60664-1:11.02

Excess voltage category: II

Permissible ambient temperature: 5 - 40 °C

Permissible relative humidity up to 31 °C: 80 %

Permissible relative humidity up to 40 °C: 50 %

Dew not permitted.

Information for the use as a medical product

Denotation: ultrasonic cleaning device

UMDNS-nomenclature (ECRI): 14-263

Purpose determination: In-vitro diagnostic sample preparation or lab process

Classified according to the directive
for in-vitro diagnostics 98/79/EG: device according to product category 5

Type, model, serial number,
year of manufacture: specifications see rating plate

Information according to IEC 601-1:

Protection class: class I

7 Disfunctions

- Oscillations are too weak and non-uniform, the oscillation noise is too loud; the cleaning effect is insufficient:
 - Has the liquid been properly degassed?
 - ⇒ sonicate for 15 minutes.
 - ⇒ change the fluid level.
 - Non-uniform noise (wobbling) is not a defect
- Minor erosion marks at the tank bottom?
 - ⇒ common wear and tear effect. The unit is faultless.
- The mains fuse can be released.
 - ⇒ The unit is equipped with an overvoltage arrester which protects it from voltage peaks; the mains fuse can be released during overcharge.
- Repairs should only be carried out by authorised service.
- Replace defective parts only by original **SONOREX** parts.

Please inform us about disfunctions in writing.

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