

Large Animal Anaesthesia Machine LAVC 2000



- ✓ with bird respirator for controlled ventilation
- ✓ Flowmeter for O₂ or N₂O (0-10 l)
- ✓ flush valve
- ✓ pressure gauge
- ✓ 10 l CO₂ absorber
- ✓ two 160 cm corrugated breathing tubes
- ✓ inhaler Y-piece with 4 adapters:(18mm / 20mm / 25mm/30mm)
- ✓ circle system for spontaneous respiration
- ✓ with 30 l breathing bag
- ✓ optional: foal set

ORDER INFORMATION:

- Code: 21 32 00 Large animal anaesthesia Machine LAVC 2000 with reconditioned Isoflurane vaporizer
- Code: 21 32 01 Large animal anaesthesia Machine LAVC 2000 with new Isoflurane vaporizer
- Code: 21 32 17 Large animal anaesthesia Machine LAVC 2000 with reconditioned Halothane vaporizer